



Seminascara filtrante FFP2M Conforme agli standard sempre DIN EN 149:2001 A1:2009 Regolamento CE 2016/425 FFP2M Filtering half mask Complies with European standards PPE EN 149:2001 A1:2009 Use directive 2016/425 *

High quality
 Made in
 Turchia

1
 PZ

IT-1. Tenere la mascherina in modo che lo stringinasso pieghevole sia posizionato nella parte superiore. 2. Posizionare i piani elastici. 3. Regolare la clip stringinasso per adattare la forma della mascherina al viso 4. Verificare la tenuta della mascherina prima di entrare in una zona pericolosa. 5. Toccare la mascherina solo sui laccetti durante la rimozione. L'esterno della mascherina può essere contaminato. Gettare immediatamente in un contenitore chiuso 8. Disinfettare le mani prima e dopo.



Mascherina FFP2 in grado di bloccare le particelle solide e i droplets (micro goccioline respiratorie) che possono trasportare virus e batteri, proteggendo le vie respiratorie. L'utilizzo di materiali d'avanguardia ha portato alla realizzazione di una mascherina efficiente, comoda e leggera. Il filtro antiparticolato offre protezione da livelli medi di polveri sottili e nebbie a base acqua ed oleosa, per il libero uso della Popolazione, per la protezione propria altrui. Il morbido bordo di tenuta facciale interno, la barra nasale adattabile a gli elastici laterali garantiscono una vestibilità confortevole per visi di conformazioni diverse. Alta filtrazione, potere filtrante tipo FFP2 > 95% -Manutenzione: Non occorre essendo monouso -Destinatari: Le maschere filtranti FFP2 devono essere utilizzate a protezione personale e per tutelare chi ci circonda

Manufacturer: AC MEDICAL
 Address: Gekturk Merkez Mah. Gekturk Cad.
 No: 3/8 Eyepozlak/ISTANBUL

IMPORTATO DA:
 ITALIAMEDICA MONDO - R&S
 mail: info@italamedica.it
 www.italamedica.it

Sul sito www.monlab.com è disponibile la dichiarazione di conformità e la scheda informativa del prodotto



2841



NR2841

Data MOD. SCADENZA
 21.01.2022 21.01.2025
 LOTTO 202221



ITALIAMEDICA
 MEDICAL PERSONAL PROTECTIVE EQUIPMENT

Turchia

AB Tip İnceleme Sertifikası EU Type-Examination Certificate

Belge No / Certificate No : 167-21-01-R01
Belgeleme Tarihi - Bir Sonraki Belge Tarihi /
Certification Date / Certificate Validity Date : 21.04.2021-23.02.2026
Belge Geçerlilik Tarihi / Document Validity Period: 5 yıl / 5 years
Firma Unvanı ve Adresi /
Company Name and Address : AC MEDİKAL
Gökürk Merkez Mah. Gökçek Cad. No:
3/A Eşişultan/ İSTANBUL

Ürün Adı / Modeller / Product Name / Models : A.C.M.I.100-02
Direktif / Directive : 2016/425 REGULATION
Modül/Kategori / Module / Category : B MODÜLÜ/ KATEGORİ III
MODULE B / CATEGORY III
Test Rapor No'ları / Test Report No : M-2021-00172, M-2021-00600

Ürün Tipi / Product Type:
- EN 149:2001+ A1:2009 Solunumla ilgili koruyucu cihazlar - Parçaçıklara karşı koruma amaçlı filtrelili yarı maskeler/ Respiratory protective devices - Filtering half masks to protect against particles

Ürünün Malzeme Bilgisi / Product Material Information: A.C.M.I.100-02 model ürünleri kumaş, elastik kayış, burun klipsi, filtre katmanı kullanılarak üretilmiştir./ A.C.M.I.100-02 model products are manufactured using fabric, elastic strap, nose clip, filter layer.

Revizyon nedeni / Reason for revision: Farklı renkte ürünler eklenmiş ve model adı revize edilmiştir./ Different colored products have been added and the model name has been revised.

Volkan AKIN
21.04.2021
Karar Verici / Approver



Öhan AKEL
21.04.2021
Şirket Müdürü / General manager









To certify the PPE product at Category III level, C2 or D module is accompanied by applying one of the conformity assessment methods along with the EU Type Examination (Module B).

Model : A.C.M.I.100-02

PPE SPECIFICATION	PERFORMANCE LEVELS
Classification	FFP2
Reusable / Single Shift Use	NR

PPE produced as a single unit to fit an individual user, all the necessary instructions for manufacturing such PPE on the basis of the approved basic model:

MARKING
MANUFACTURER: AC MEDİKAL
PPE TYPE : - EN 149:2001+ A1:2009 Respiratory protective devices - Filtering half masks to protect against particles
MODEL: A.C.M.I.100-02
PICTOGRAM AND PERFORMANCE LEVELS: EN 149:2001+ A1:2009 FFP2 NR       NB 2841 Or Condition of Storage

MNA LABORATORIES SAN. TIC. LTD. ŞTİ declares that the above-mentioned product meets the requirements of the directive according to the EU Directive 2016/425, the safety of the product is covered by the conditions and use specified in this certificate and in the technical file.

ATTACHMENTS (167-21-01-R01)



DOCUMENTS IN THE TECHNICAL FILE

- Basic Health Safety Requirements
- Risk Assessment
- Test Reports
- Technical Report

Report No : 167-21-01-R01

Report Date : 21.04.2021

Application No : 167-21-01

- COMPANY INFORMATION:**
AC MEDICAL
Göktürk Merkez Mah. Göktürk Cad. No: 3/A Eyüpsultan/ İSTANBUL
Tel: +90 532 415 67 27
- PPE INFORMATION:**
Disposable and non-sterile half mask made of particulate protection filter material.
- PPE TYPE IDENTIFICATION**
EN 149:2001+A1:2009 Respiratory protective devices – Filtering half masks to protect against particles - Requirements, testing, marking
- PPE PICTURES**



5. PPE DIMENSIONS:

A.C.M.I.100-02 model has been found to be produced using standart sizes.

6. PPE PRODUCT MATERIAL INFORMATION:

The product is made of elastic strap, nonwoven fabric on the outer and inner layers and filter material on the middle layer.

7. ESSENTIAL HEALTH AND SAFETY REQUIREMENTS

- A visual inspection was made according to EN 149:2001 +A1:2009 for ergonomics.
- Protection levels and degrees are defined by the manufacturer.
- Suitable construction materials were determined by visual inspection according to EN 149:2001 +A1:2009.

8. ANALYSIS AND EVALUATIONS:

EN 149:2001 +A1:2009

TESTS	PARAMETER	PERFORMANCE LEVELS			RESULTS	PERFORMANCE LEVELS	EVALUATION
		FFP1	FFP2	FFP3			
Banned Azo Dyes	< 30 mg/kg				< 5 mg/kg	-	PASS
Part 7.3 Visual inspection	Shall also the marking and the information supplied by the manufacturer				Appropriate	-	PASS
Part 7.4 Packaging	Particle filtering half mask shall be offered for sale packaged in such a way that they are protected against mechanical damage and contamination before use.				Appropriate	-	PASS
Part 7.5 Material	When conditioned in accordance 8.3.1 & 8.3.2 the particle filter half mask shall not collapse.				Appropriate	-	PASS
Part 7.6 Cleaning and disinfecting	After cleaning and disinfecting the re-usable particle filtering half mask shall satisfy the penetration requirement of the relevant class.				Not applicable	-	Not applicable
Part 7.7 Practical performance	No negative comments should be made by the test subject regarding any of the criteria evaluated.				Appropriate	-	PASS
Part 7.8 Finish of parts	Parts of the device likely to come into contact with the wearer shall have no sharp edge or burrs.				Appropriate	-	PASS
TESTS	PARAMETER	PERFORMANCE LEVELS			RESULTS	PERFORMANCE LEVELS	EVALUATION
		FFP1	FFP2	FFP3			
Part 7.9.1 Total inward leakage	At least 46 out of the 50 individual exercise result	<25	<11	<5	See the table below	FFP2	PASS
	At least 8 out of the 10 individual wearer arithmetic means	<22	<8	<2	See the table below	FFP2	PASS

	Total Inward Leakage (%)					
	Exercise 1	Exercise 2	Exercise 3	Exercise 4	Exercise 5	Average
Subject 1 (As received)	6.2	7.1	5.3	6.6	6.7	6.4
Subject 2 (As received)	8.0	5.2	6.2	8.2	6.5	6.8
Subject 3 (As received)	7.7	5.3	5.8	6.5	6.4	6.3
Subject 4 (As received)	7.4	5.2	5.9	8.2	7.7	6.9
Subject 5 (As received)	7.3	8.0	7.8	8.3	9.0	8.1
Subject 6 (After temperature conditioning)	7.1	8.4	7.7	7.3	7.2	7.5
Subject 7 (After temperature conditioning)	7.4	5.3	5.9	6.5	8.7	6.8
Subject 8 (After temperature conditioning)	6.3	7.6	7.3	6.3	7.2	6.9
Subject 9 (After temperature conditioning)	6.3	8.3	7.1	7.2	7.4	7.3
Subject 10 (After temperature conditioning)	6.1	7.0	8.5	8.2	8.4	7.6

Subject facial dimensions

Subject	Face Length (mm)	Face Width (mm)	Face Depth (mm)	Mouth Width (mm)
1	133	132	132	65
2	125	144	116	67
3	126	135	124	75
4	123	133	134	74
5	117	135	122	73
6	122	142	133	66
7	113	132	114	75
8	135	123	123	65
9	122	135	133	74
10	135	142	125	83

TESTS	PARAMETER	PERFORMANCE LEVELS			RESULTS	PERFORMANCE LEVELS	EVALUATION
		FFP1	FFP2	FFP3			
Part 7.9.2 Penetration of filter material	Sodium chloride, 95 L/min % max	% 20	% 6	% 1	See the table below	FFP2	PASS
	Paraffin oil, 95 L/min % max	% 20	% 6	% 1	See the table below	FFP2	PASS

Penetration of filter material	Sodium Chloride (%)	Paraffin Oil (%)
As received	3.5	3.6
As received	3.4	4.2
As received	4.1	3.7
After the simulated wearing treatment	3.6	4.1
After the simulated wearing treatment	4.2	3.9
After the simulated wearing treatment	3.6	3.7
Mechanical strength and temperature conditioning	4.8	5.1
Mechanical strength and temperature conditioning	4.4	4.7
Mechanical strength and temperature conditioning	4.6	4.9

TESTS	PARAMETER	PERFORMANCE LEVELS			RESULTS	PERFORMANCE LEVELS	EVALUATION
		FFP1	FFP2	FFP3			
Part 7.10 Compatibility with skin	Materials shall not be known to be likely to cause irritation or any other adverse effect to health				Appropriate	-	PASS
Part 7.11 Flammibility	Mask shall not burn or not to continue to burn for more than 5 s				Flame not seen	-	PASS
Part 7.12 Carbondioxide content of the inhalation air	Shall not exceed an average of % 1				0,87 0,82 0,85	-	PASS
Part 7.13 Head harness	It can be donned and removed easily				Appropriate	-	PASS
Part 7.14 Field of vision	The field of vision shall acceptable in practical performance test.				Appropriate	-	PASS
Part 7.15 Exhalation valve(s)	It shall withstand axially a tensile force of 10 N apply for 10 s. If fitted, shall continue to operate correctly after a continuous exhalation flow of 300 L/min over a period of 30 s.				Not applicable	-	Not applicable

TESTS	PARAMETER	PERFORMANCE LEVELS			RESULTS	PERFORMANCE LEVELS	EVALUATION
		FFP1	FFP2	FFP3			
Part 7.16 Breathing Resistance	Inhalation 30L/min	0,6 mbar	0,7 mbar	1,0 mbar	See the table below	FFP2	PASS
	Inhalation 95L/min	2,1 mbar	2,4 mbar	3,0 mbar	See the table below	FFP2	PASS
	Exhalation 160L/min	3,0 mbar	3,0 mbar	3,0 mbar	See the table below	FFP2	PASS

Breathing Resistance (mbar)	Inhalation 30L/min	Inhalation 95L/min
As received	0,4	1,7
As received	0,5	1,8
As received	0,4	1,8
After temperature conditioning	0,5	1,7
After temperature conditioning	0,5	1,8
After temperature conditioning	0,5	1,7
After the simulated wearing treatment	0,5	1,7
After the simulated wearing treatment	0,4	1,8
After the simulated wearing treatment	0,4	1,8

Breathing Resistance 160L/min (mbar)	Facing directly ahead	Facing vertically upwards	Facing vertically downwards	Lying on the left side	Lying on the right side
As received	2,3	2,3	2,2	2,3	2,3
As received	2,2	2,3	2,3	2,3	2,3
As received	2,3	2,2	2,3	2,3	2,3

After temperature conditioning	2,2	2,3	2,3	2,2	2,3
After temperature conditioning	2,2	2,2	2,3	2,3	2,3
After temperature conditioning	2,3	2,2	2,3	2,2	2,3
After the simulated wearing treatment	2,3	2,3	2,3	2,3	2,3
After the simulated wearing treatment	2,3	2,3	2,3	2,2	2,3
After the simulated wearing treatment	2,3	2,3	2,3	2,3	2,3

TESTS	PARAMETER	PERFORMANCE LEVELS			RESULTS	PERFORMANCE LEVELS	EVALUATION
		FFP1	FFP2	FFP3			
Part 7.17 Clogging	After clogging the inhalation resistances shall not exceed. (valved)	4 mbar	5 mbar	7 mbar	Not applicable	-	Not applicable
	The exhalation resistance shall not exceed 3 mbar at 160 L/ min continuous flow. (valved)				Not applicable	-	Not applicable
	After clogging the inhalation and exhalation resistances shall not exceed. (valveless)	3 mbar	4 mbar	5 mbar	Not applicable	-	Not applicable
Part 7.18 Demountable part	All demountable parts (if fitted) shall be readily connected and secured were possible by hand.				Not applicable	-	Not applicable

9. DECISION PROPOSAL

Analysis and examinations A.C M.L100-02 model coded personal protective equipment; Respiratory Protective Devices EN 149:2001 +A1:2009- Filtered Half Masks for Protection Against Particles - Properties, Experiments and Marking standards are evaluated. It is recommended to be certified at the performance levels specified as a result of technical evaluations.

10. ATTACHMENTS

- Basic Health Safety Requirements
- Risk Assessment
- User Instruction

Reason for revision : Different colored products have been added and the model name has been revised.

CONTROLLER : VOLKAN AKIN

SIGN :

DATE : 21.04.2021



Notified Body Number: 2841

Report No : 167-21-01-R01-01

Report Date : 21.04.2021

Application No : 167-21-01-R01-01

1. COMPANY INFORMATION:

AC MEDİKAL

Göktürk Merkez Mah. Göktürk Cad. No: 3/A Eyüpsultan/ İSTANBUL

Tel: +90 532 415 67 27

2. PPE INFORMATION:

Disposable and non-sterile half mask made of particulate protection filter material.

3. PPE TYPE IDENTIFICATION

EN 149:2001+A1:2009 Respiratory protective devices – Filtering half masks to protect against particles - Requirements, testing, marking

4. PPE PICTURES



A.C M.I.100-02

U-FRM-056-REV.00.YAYIN TARİHİ:20.11.2020

Notified Body Number: 2841

5. PPE DIMENSIONS:

A.C M.I.100-02 model has been found to be produced using standard sizes.

6. PPE PRODUCT MATERIAL INFORMATION:

The mask is made of elastic strap, nonwoven fabric on the outer and inner layers and filter material on the middle layer.

7. ESSENTIAL HEALTH AND SAFETY REQUIREMENTS

- A visual inspection was made according to EN 149:2001 +A1:2009 for ergonomics.
- Protection levels and degrees are defined by the manufacturer.
- Suitable construction materials were determined by visual inspection according to EN 149:2001 +A1:2009.

8. ANALYSIS AND EVALUATIONS:

EN 149:2001 +A1:2009

TESTS	PARAMETER	PERFORMANCE LEVELS			RESULTS	PERFORMANCE LEVELS	EVALUATION
		FFP1	FFP2	FFP3			
Banned Dyes	Azo	< 30 mg/kg			< 5 mg/kg	-	PASS
Part 7.3 Visual inspection	Shall also the marking and the information supplied by the manufacturer			Appropriate	-	PASS	
Part 7.4 Packaging	Particle filtering half mask shall be offered for sale packaged in such a way that they are protected against mechanical damage and contamination before use.			Appropriate	-	PASS	
Part 7.5 Material	When conditioned in accordance 8.3.1 & 8.3.2 the particle filter half mask shall not collapse.			Appropriate	-	PASS	
Part 7.6 Cleaning and disinfecting	After cleaning and disinfecting the re-usable particle filtering half mask shall satisfy the penetration requirement of the relevant class.			Not applicable	-	Not applicable	
Part 7.7 Practical performance	No negative comments should be made by the test subject regarding any of the criteria evaluated.			Appropriate	-	PASS	
Part 7.8 Finish of parts	Parts of the device likely to come into contact with the wearer shall have no sharp edge or burrs.			Appropriate	-	PASS	
TESTS	PARAMETER	PERFORMANCE LEVELS			RESULTS	PERFORMANCE LEVELS	EVALUATION
		FFP1	FFP2	FFP3			
Part 7.9.1 Total inward leakage	At least 46 out of the 50 individual exercise result	<25	<11	<5	See the table below	FFP2	PASS
		<22	<8	<2			

U-FRM-056-REV.00.YAYIN TARİHİ:20.11.2020

Total Inward Leakage (%)						
	Exercise 1	Exercise 2	Exercise 3	Exercise 4	Exercise 5	Average
Subject 1 (As received)	6.5	6.2	5.8	8.1	8.2	7.0
Subject 2 (As received)	6.0	6.4	5.2	8.3	8.1	6.8
Subject 3 (As received)	7.3	6.3	7.9	7.5	7.9	7.4
Subject 4 (As received)	7.0	8.0	7.8	9.2	8.1	8.0
Subject 5 (As received)	6.7	8.2	8.5	4.1	8.1	7.1
Subject 6 (After temperature conditioning)	6.6	7.3	7.1	6.6	8.1	7.1
Subject 7 (After temperature conditioning)	6.4	7.6	7.0	7.9	6.5	7.1
Subject 8 (After temperature conditioning)	6.7	7.0	6.7	8.1	8.0	7.3
Subject 9 (After temperature conditioning)	6.7	6.9	6.6	7.9	6.5	6.9
Subject 10 (After temperature conditioning)	6.6	7.7	7.9	7.7	6.5	7.3

Subject facial dimensions

Subject	Face Length (mm)	Face Width (mm)	Face Depth (mm)	Mouth Width (mm)
1	133	132	132	65
2	125	144	116	67
3	126	135	124	75
4	123	133	134	74
5	117	135	122	73
6	122	142	133	66
7	113	132	114	75
8	135	123	123	65
9	122	135	133	74
10	135	142	125	83

TESTS	PARAMETER	PERFORMANCE LEVELS			RESULTS	PERFORMANCE LEVELS	EVALUATION
		FFP1	FFP2	FFP3			
Part 7.9.2 Penetration of filter material	Sodium chloride, 95 L/min % max	% 20	% 6	% 1	See the table below	FFP2	PASS
	Paraffin oil, 95 L/min % max	% 20	% 6	% 1	See the table below	FFP2	PASS
Penetration of filter material					Sodium Chloride (%)	Paraffin Oil (%)	
	As received				3.6	3.9	
	As received				3.5	4.5	
	As received				4.1	4.1	
	After the simulated wearing treatment				3.6	4.0	
	After the simulated wearing treatment				3.3	4.2	
	After the simulated wearing treatment				3.0	4.6	
	Mechanical strength and temperature conditioning				4.9	5.0	
	Mechanical strength and temperature conditioning				5.1	5.1	
	Mechanical strength and temperature conditioning				5.0	5.4	

TESTS	PARAMETER	PERFORMANCE LEVELS			RESULTS	PERFORMANCE LEVELS	EVALUATION
		FFP1	FFP2	FFP3			
Part 7.10 Compatibility with skin	Materials shall not be known to be likely to cause irritation or any other adverse effect to health				Appropriate	-	PASS
Part 7.11 Flammability	Mask shall not burn or not to continue to burn for more than 5 s				Flame not seen	-	PASS
Part 7.12 Carbon dioxide content of the inhalation air	Shall not exceed an average of % 1				0,88 0,83 0,89	-	PASS
Part 7.13 Head harness	It can be donned and removed easily				Appropriate	-	PASS
Part 7.14 Field of vision	The field of vision shall acceptable in practical performance test.				Appropriate	-	PASS
Part 7.15 Exhalation valve(s)	It shall withstand axially a tensile force of 10 N apply for 10 s. If fitted, shall continue to operate correctly after a continuous exhalation flow of 300 L/min over a period of 30 s.				Not applicable	-	Not applicable

TESTS	PARAMETER	PERFORMANCE LEVELS			RESULTS	PERFORMANCE LEVELS	EVALUATION
		FFP1	FFP2	FFP3			
Part 7.16 Breathing Resistance	Inhalation 30L/min	0,6 mbar	0,7 mbar	1,0 mbar	See the table below	FFP2	PASS
	Inhalation 95L/min	2,1 mbar	2,4 mbar	3,0 mbar	See the table below	FFP2	PASS
	Exhalation 160L/min	3,0 mbar	3,0 mbar	3,0 mbar	See the table below	FFP2	PASS

Breathing Resistance (mbar)		Inhalation 30L/min	Inhalation 95L/min
	As received	0.4	2.2
	As received	0.5	2.2
	As received	0.5	2.1
	After temperature conditioning	0.4	2.1
	After temperature conditioning	0.5	2.1
	After temperature conditioning	0.4	2.2
	After the simulated wearing treatment	0.4	2.2
	After the simulated wearing treatment	0.5	2.1
	After the simulated wearing treatment	0.4	2.1

Breathing Resistance 160L/min (mbar)	Facing directly ahead	Facing vertically upwards	Facing vertically downwards	Lying on the left side	Lying on the right side
	As received	2.8	2.6	2.8	2.8
As received	2.8	2.7	2.8	2.8	2.7
As received	2.8	2.6	2.7	2.8	2.7

After temperature conditioning	2,7	2,8	2,7	2,8	2,8
After temperature conditioning	2,8	2,8	2,8	2,8	2,8
After temperature conditioning	2,8	2,8	2,8	2,7	2,8
After the simulated wearing treatment	2,8	2,8	2,8	2,8	2,7
After the simulated wearing treatment	2,7	2,8	2,7	2,8	2,8
After the simulated wearing treatment	2,8	2,8	2,8	2,7	2,8

TESTS	PARAMETER	PERFORMANCE LEVELS			RESULTS	PERFORMANCE LEVELS	EVALUATION
		FFP1	FFP2	FFP3			
Part 7.17 Clogging	After clogging the inhalation resistances shall not exceed. (valved)	4 mbar	5 mbar	7 mbar	Not applicable	-	Not applicable
	The exhalation resistance shall not exceed 3 mbar at 160 L/ min continuous flow. (valved)				Not applicable	-	Not applicable
	After clogging the inhalation and exhalation resistances shall not exceed. (valveless)	3 mbar	4 mbar	5 mbar	Not applicable	-	Not applicable
Part 7.18 Demountable part	All demountable parts (if fitted) shall be readily connected and secured were possible by hand.				Not applicable	-	Not applicable

9. DECISION

Analysis and examinations A.C. M.L100-02 model coded personal protective equipment; Respiratory Protective Devices EN 149:2001 +A1:2009- Filtered Half Masks for Protection Against Particles - Properties, Experiments and Marking standards are evaluated. The homogeneity of the production was monitored at the performance levels determined as a result of the technical evaluations made within the scope of MODULE C2.

10. ATTACHMENTS

- Basic Health Safety Requirements
- Risk Assessment
- Test Reports (M-2021-00600)
- User Instruction

CONTROLLER : VOLKAN AKIN
SING :
DATE : 21.04.2021



EU DECLARATION OF CONFORMITY

MANUFACTURER
AC MEDİKAL

Göktürk Merkez Mah. Göktürk Cad. No:3/A Eyişultam' İSTANBUL / TURKEY

PRODUCT DESCRIPTION

Brand Name: AC MASK Model: AC M.L 100-02
Particle Filtering Half Mask
Class: FFP2 NR

Particle Filtering Half Face Mask in Category III product according to (EU) 2016/425 Personal Protective Equipment Regulation

The Manufacturer declares on his sole responsibility that the product above is, under conditions of normal use and conditions defined by the Manufacturer, safe and meets all the necessary legal conditions and requirements. The product is a personal protective equipment that is intended for single use and solely in accordance with the Manufacturer's instructions.

The Conformity is ensured with the following mechanism:

- Complies with EU 2016/425 Personal Protective Equipment Regulation establishing technical requirements for Category III products,
- Complies with Essential Health and Safety Requirements of Technical harmonised standard EN 149:2001 +A1:2009
- All required tests referred in above standards are conducted,
- Complies with other relevant harmonized legislation and community standards
- For the assessment of conformity the EU Type Examination certificate (Serial No: 167-21-01-R01) is issued, after all technical evaluations for conformity to the regulation and harmonised standards conducted, by:
 - MNA Laboratuvarları San. Tic. Ltd. Şti., as Notified Body number 2841
- The product is under surveillance of some Notified Body, NB 2841 according to the Annex III (Module C2) of the PPE Regulation (EU) 2016/425, for quality assurance.

MARKING, LABELLING

Marking, labelling and user information are prepared in accordance with EU 2016/425 Personal Protective Equipment Regulation and the harmonised product standards given above.

MEASURES TO ENSURE CONFORMITY

The Manufacturer declares that he has taken all necessary measures to ensure the conformity of products placed on the market with technical documentation and technical requirements for this type of product.

AYKUT CENBERCİ
General Manager
22/12/2021

AC MEDİKAL
AYKUT CENBERCİ
Göktürk Merkez Mah. Göktürk Cad. No:3/A
Eyişultam' İstanbul / Türkiye
Güvenli İşleri / 0212 398 70 02 28



2841